

turns more than heads



Premium League



Zentrifugen
Centrifuges

Produktfamilie | Product family



Sigma 1-14



–

0,2 - 2 ml

48 ml

14.800

16.163



Sigma 1-14K

0,2 - 2 ml

0,2 - 2 ml

48 ml

15.000

16.602



Sigma 1-16



–

0,2 - 5 ml

72 ml

15.000

20.627



Sigma 1-16K




–

0,2 - 5 ml

72 ml

15.000

20.627



Sigma 3-30KS

0,2 - 100 ml; MTP

0,2 - 94 ml

564 ml

30.000

70.121



Sigma 4-16S

0,2 - 750 ml; MTP

0,2 - 250 ml

3.000 ml

13.500

21.917



Sigma 4-16KS

0,2 - 750 ml; MTP

0,2 - 250 ml

3.000 ml

15.000

25.155



Sigma 4-5KRL

1,1 - 15 ml; MTP

–

1.500 ml

4.700

4.643



Symbol Typ / Type



Ausschwingrotor
Swing-out rotor



Winkelrotor
Fixed-angle rotor



max. Kapazität
max. capacity



Sigma 1-7

1,1 - 15 ml

–

90 ml

8.000

6.153



Sigma 2-7

0,2 - 100 ml; MTP, Cyto

5 - 15 ml

450 ml

4.000

2.540



Sigma 2-16P

0,2 - 100 ml; MTP

0,2 - 50 ml

450 ml

15.000

20.627



Sigma 2-16KL


0,2 - 100 ml; MTP

0,2 - 85 ml

450 ml

15.300

21.913



Sigma 4-5L

0,2 - 750 ml; MTP

–

3.000 ml

4.700

4.643



Sigma 4-5KL

0,2 - 750 ml; MTP

2 - 94 ml

3.000 ml

13.500

20.376



Sigma 6-16S

0,2 - 1.000 ml; MTP, BB

0,2 - 500 ml

4.000 ml

13.500

21.917



Sigma 6-16KS

0,2 - 1.000 ml; MTP, BB

0,2 - 500 ml

4.000 ml

15.000

25.419



max. Drehzahl, min^{-1}
max. speed, rpm



max. RZB, x g
max. RCF, x g



Kühlzentrifuge
Refrigerated centrifuge



Medizinprodukt für IVD Anwendungen
Medical device for IVD applications



Sigma 3-16L

0,2 - 400 ml; MTP

0,2 - 94 ml

1.600 ml

15.300

21.460

IVD



Sigma 3-16KL

0,2 - 400 ml; MTP

0,2 - 94 ml

1.600 ml

15.300

21.913

IVD



Sigma 3-18KS

0,2 - 400 ml; MTP

0,2 - 94 ml

1.600 ml

18.000

30.070

IVD



Sigma 8KS

0,2 - 1.500 ml; MTP

0,2 - 1.000 ml

12.000 ml

10.500

20.461

IVD

Sigma 8KS Clinic

BBS, 1.500 ml, MTP

12 Blood bag systems

4.100

5.394

IVD

Sigma 8KBS

12 Blutbeutelssysteme
Blood bag systems

12 Blood bag systems

5.394

MD

More
Information
Please find details and
brochures on our website:





Spincontrol Basic



Spincontrol L

Spincontrol L
Mehr Infos erhalten
Sie in unserem Video:

Please watch our video:



Spincontrol S

Spincontrol S
Mehr Infos erhalten
Sie in unserem Video:

Please watch our video:





Sigma
Laborzentrifugen GmbH
An der Unteren Söse 50
37520 Osterode • Germany
Tel. +49 (0) 55 22 50 07-0
Fax +49 (0) 55 22 50 07-12
info@sigma-zentrifugen.de
www.sigma-zentrifugen.de