



### **Sigma 4-5L – High Capacity meets Power and Silence**

High throughput performance is the domain of the Sigma 4-5L laboratory centrifuge, which is dedicated to the use with a high capacity swing-out rotor. The new drive system has been optimized specifically in three key areas:

- high motor torque enabling fast acceleration
- minimum temperature increase
- minimum noise during operation

The result is an outstandingly powerful yet quiet instrument. Even more comfort for operating the 4-5L is added by the easy to clean stainless steel bowl, the automated lid opening and the single-knob control.

#### **Aerodynamically shaped air intake**



*Air flow pattern simulation has been applied for achieving optimum effectiveness of the air circulation through the instrument.*

#### **SIGMA Laborzentrifugen GmbH**

An der Unteren Söse 50  
37520 Osterode  
Germany

Tel.: +49 5522 5007-0  
Fax: +49 5522 5007-12  
info@sigma-zentrifugen.de  
www.sigma-zentrifugen.de