

New "Cleanroom" Accessory Range from SIGMA

Robust, Chemical Resistant, Value-Preserving

Long durability and high chemical resistance of aluminium accessories are essential for safety at work.

For the valuable aluminium buckets and rotors SIGMA is launching a new, extremely resistant coating. With this new coating the accessories are protected against corrosion – especially when working with corrosive chemicals. Even caustic soda and hydrochloric acids have no chance against the "Cleanroom" coating from SIGMA.

The aluminium accessories with the new special coating have a high durability. The surfaces are also food-safe (FDA certified) and easy to clean without leaving any residue.

Keeping it Clean

The SIGMA "Cleanroom" surface is especially engineered for clean rooms.

Due to the high resistance of the coating there are no interfering particles produced by abrasion during centrifugation. The robust coating prevents the occurrence of particle wear debris thus eliminating associated product contamination.



SIGMA Laborzentrifugen GmbH

An der Unteren Söse 50
37520 Osterode
Germany

Tel.: +49 5522 5007-0
Fax: +49 5522 5007-12
info@sigma-zentrifugen.de
www.sigma-zentrifugen.de