

Premium League



Zentrifugen
Centrifuges

Produktfamilie / Product family



Sigma 1-14




0,2 - 2 ml

48 ml

14.800

16.163

Sigma 1-14K



0,2 - 2 ml

0,2 - 2 ml

48 ml

15.000

16.602



Sigma 1-16



0,2 - 5 ml

72 ml

15.000

20.627



Sigma 1-16K




0,2 - 5 ml

72 ml

15.000

20.627



Sigma 3-30KS

0,2 - 100 ml; MTP

0,2 - 94 ml

564 ml

30.000

70.121



Sigma 4-16S

0,2 - 750 ml; MTP

0,2 - 250 ml

3.000 ml

13.500

20.376



Sigma 4-16KS

0,2 - 750 ml; MTP

0,2 - 250 ml

3.000 ml

15.000

25.155





Sigma 1-7

0,2 - 15 ml

–

90 ml

8.000

6.153



Sigma 2-7

0,2 - 100 ml; MTP, Cyto

5 - 15 ml

450 ml

4.000

2.540



Sigma 2-16P

0,2 - 100 ml; MTP

0,2 - 50 ml

450 ml

15.000

20.627



Sigma 2-16KL

0,2 - 100 ml; MTP

0,2 - 85 ml

450 ml

15.300

21.913



500



500



650



800




1000




1000



1300



1500



MTP



BB



Sigma 4-5L

0,2 - 750 ml; MTP

–

3.000 ml

4.700

4.643



Sigma 4-5KRL

1,1 - 15 ml; MTP

–

1.500 ml

4.700

4.470



Sigma 6-16S

0,2 - 1.000 ml; MTP, BB

0,2 - 500 ml

4.000 ml

13.500

20.376



Sigma 6-16KS

0,2 - 1.000 ml; MTP, BB

0,2 - 500 ml

4.000 ml

15.000

25.419





Sigma 3-16L

0,2 - 400 ml; MTP

0,2 - 94 ml

1.600 ml

15.300

21.460

IVD



Sigma 3-16KL

0,2 - 400 ml; MTP

0,2 - 94 ml

1.600 ml

15.300

21.913

IVD



Sigma 3-18KS

0,2 - 400 ml; MTP

0,2 - 94 ml

1.600 ml

18.000

30.070

Refrigerated



Sigma 8KS / 8KBS

0,2 - 1.500 ml; MTP, BB

0,2 - 1.000 ml


12.000 ml

10.500

20.461

Refrigerated

Blut- und Uringefäße
Blood collection tubes



Spitzgefäße
Conical tubes



Rundbodengefäße
Round bottom tubes



Symbol Typ / Type

Ausschwingrotor
Swing-out rotor

Winkelrotor
Fixed-angle rotor

max. Kapazität
max. capacity

max. Drehzahl, min⁻¹
max. speed, rpm

max. RZB, x g
max. RCF, x g

Kühlzentrifuge
Refrigerated centrifuge

IVD Medizinprodukt für IVD Anwendungen
Medical device for IVD applications

Spincontrol family



Spincontrol Basic



Spincontrol L



Spincontrol S





Sigma
Laborzentrifugen GmbH
An der Unteren Söse 50
37520 Osterode • Germany
Tel. +49 (0) 55 22 50 07-0
Fax +49 (0) 55 22 50 07-12
info@sigma-zentrifugen.de
www.sigma-zentrifugen.de